


Minimalism in Sinus

Easy & Safe

Shorter
Chairtime

DASK *Simple*








Minimalism in Sinus

DASK *Simple*


Dentium

Bone Compaction

Safe membrane elevation



50rpm 70N·cm
Non irrigation




Detaching



Sinus Lifting



Round Instrument



SPECIFICATION



Compaction Drill



Top view



Diameter
Ø 3.0

※Drills be replaced after 50 uses

Diameter	L	Art. No
Ø 3.0	35.0mm	SRT3035AS
Ø 3.5	35.0mm	SRT3535AS


Lateral Drill



※Drills be replaced after 40 uses

Diameter	L	Art. No
Ø 4.0	28.0mm	SLD0428
Ø 4.0	40.0mm	SLD0440

Bone Spreader




Diameter	L	Art. No
Ø 5.0	33.5mm	SSP33

DASK *Simple*

Lateral approach

NEW



✓ Handy easy

Unlimited approach with Lateral Drill for Sinus solution

✓ Free crestal and lateral excess

1,000rpm 30~45N·cm with irrigation

