

The Trust

Implant Stability Measuring Device

- Usable at any angle
- Similar Osstell ISQ
- Disposable vinyl cap



FEATURES




Easy operation



Quick check result
(7 times in 3 seconds)



LCD window



Dentium ISV Scale

Implant Stability	Dentium ISV Scale					ISV
	LOW STABILITY		MEDIUM STABILITY		HIGH STABILITY	
Surgical protocol	1- or 2-stage	1- or 2-stage	1-stage	1-stage	1-stage	Implant Stability Value
Restorative protocol	Traditional loading	Early loading (6-8w)	Early loading(2-4w)	Immediate loading		

Dentium

The Trust

Implant Stability Measuring Device

- Usable at any angle
- Similar Osstell ISQ
- Disposable vinyl cap



FEATURES




Easy operation



Quick check result
(7 times in 3 seconds)



LCD window



Dentium ISV Scale

Implant Stability	Dentium ISV Scale					ISV
	LOW STABILITY		MEDIUM STABILITY		HIGH STABILITY	
Surgical protocol	1- or 2-stage	1- or 2-stage	1-stage	1-stage	1-stage	Implant Stability Value
Restorative protocol	Traditional loading	Early loading (6-8w)	Early loading(2-4w)	Immediate loading		

Dentium