

bright *Chair* 2

Brochure



Dentium

bright Chair 2



Clear Visual Field

Up to 50,000 LUX, ensuring clarity and a wide view with a broad FOV




Blue Light
Filtering Function

Maximized Spatial Efficiency


Compact size allows use of remaining space when installing chair

bright chair 2



4m

Other Companies



5m




bright chair 2
Shortcut to
Promotional Video




Worthy Long Lasting Tech

Robust floor mount type structure allows for long-term use without minor breakdowns




One Person Treatment

Specially designed self-section module for easy single user treatment




Chair Dimension


Chair Out Size





Chair Up




Chair Down




Chair Out Size



Backrest Up



Backrest Down



Specification

Lift System	Hydraulic
Lift Height	490mm ~ 750mm (250mm)
Recline Range	0° ~ 65° (65°)
Headrest	Two-Joint
Dental Light	LED Light
Controller	Switch + LCD
Spittoon	Manual Rotation

Seat Color

Most Selected



Indigo Navy



Ivory Beige



Cherry Red



Ocean Blue



Olive Green



Dark Cyan



Royal Wine



Mocha Brown



Solid Black



Indigo Navy



Dark Cyan



Ocean Blue



Mocha Brown

Instrument

Line up (Basic)

High Handpiece	Low Handpiece	Coupling
 M800L(NSK) Optic - Power: 23W - Speed: 380,000~450,000min⁻¹ - Head Size: ø 10.6 x H12.4 mm	   FX205 M4/FX65/FX25(NSK) Non-Optic Midwest 4 hole · Max: 25,000 min⁻¹	 PTL-CL-LED(NSK) Optic - For NSK: Stainless Body · Titanium Coating · LED Coupling  KCL-LED(KAVO) Optic - For KaVo® MULTIflex® LUX with water volume adjuster - Stainless Body · Titanium Coating · LED Coupling
 M900L(NSK) Optic - Power: 26W - Speed: 325,000~430,000min⁻¹ - Head Size: ø 10.8 x H12.1 mm	Scaler  bright Scaler(Dentium) 28kHz~32kHz · 71g  Varios 170(NSK)  dmetec	
 M800KL(NSK) Optic - Power: 23W - Speed: 380,000~450,000min⁻¹ - Head Size: ø 10.6 x H12.4 mm		
 M900KL(NSK) Optic - Power: 26W - Speed: 325,000~430,000min⁻¹ - Head Size: ø 12.1 x H13.3 mm		

Option

High Handpiece	Angle	Motor
 Z900L(NSK) Optic - Power: 23 W, Speed: 360,000 ~ 440,000 min⁻¹ - Head Size: ø 12.5 x H13.1mm	  M65/M25(NSK) Non-Optic M25L(NSK) Optic - Stainless Body · Backflow Prevention Valve - Midwest 4 hole · Max: 25,000 min⁻¹  M95L(NSK) Optic 4-point Douche Spray · For FG Bur (ø1.6) Ceramic Ball Bearing - Max Speed: 200,000mi	 bright Elec(Dentium) - Speed: 40,000rpm ± 10% - Torque: 3.0Ncm ± 10% - Weight: 62g - Light: Led light, 10,000lux  M205(NSK) - Stainless Body · Backflow Prevention Valve · Midwest 4 hole - Max: 25,000 min⁻¹  NLX nano(NSK) - Speed: 1,000 - 40,000 r/min (NLXBF) - Torque: Max. 3.4Ncm - Weight: 61g (micromotor only) - Light: White LED
 Z800L(NSK) Optic - Power: 23 W, Speed: 360,000 ~ 440,000 min⁻¹ - Head Size: ø 10.8 x H12.1mm		

Dentium

Dentium



This publication is for international distribution. It may differ from the certification and sales products in some countries.

BRC2_2023.08 [Rev.0]

76, Changyong-daero 256beon-gil, Yeongtong-gu, Suwon-si, Gyeonggi-do, Korea 16229
E-mail : overseas@dentium.com / www.globaldentium.com